

Keratin, type I cytoskeletal 9

Catalog No: #AP85333



Package Size: #AP85333-1 50ug #AP85333-2 100ug #AP85333-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Keratin, type I cytoskeletal 9 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Dog |
| Other Names | Cytokeratin-9,CK-9,Keratin-9,K9 |
| Accession No. | Uniprot ID:O18740 |
| Uniprot | O18740 |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:KRT9

Note: This product is for in vitro research use only