

Melanoma-associated antigen 1

Catalog No: #AP85336



Package Size: #AP85336-1 50ug #AP85336-2 100ug #AP85336-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Melanoma-associated antigen 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Antigen MZ2-E,Cancer/testis antigen 1.1,CT1.1,MAGE-1 antigen,MAGE1,MAGE1A
Accession No.	Uniprot ID:P43355
Uniprot	P43355
GeneID	4100;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MAGEA1

Note: This product is for in vitro research use only