

NADH-cytochrome b5 reductase 3

Catalog No: #AP85343



Package Size: #AP85343-1 50ug #AP85343-2 100ug #AP85343-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NADH-cytochrome b5 reductase 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	B5R,Diaphorase-1,DIA1
Accession No.	Uniprot ID:P07514
Uniprot	P07514
GeneID	515773;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CYB5R3

Note: This product is for in vitro research use only