

Kin of IRRE-like protein 1

Catalog No: #AP85344



Package Size: #AP85344-1 50ug #AP85344-2 100ug #AP85344-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Kin of IRRE-like protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Kin of irregular chiasm-like protein 1,Nephrin-like protein 1,KIRREL1,NEPH1
Accession No.	Uniprot ID:Q96J84
Uniprot	Q96J84
GeneID	55243;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:KIRREL

Note: This product is for in vitro research use only