

Lymphatic vessel endothelial hyaluronic acid receptor 1

Catalog No: #AP85346

Orders: order@signalwayantibody.com

Package Size: #AP85346-1 50ug #AP85346-2 100ug #AP85346-3 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Lymphatic vessel endothelial hyaluronic acid receptor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	LYVE-1,Cell surface retention sequence-binding protein 1,CRSBP-1,Extracellular link domain-containing protein 1,CRSBP1,XLKD1
Accession No.	Uniprot ID:Q6UC88
Uniprot	Q6UC88
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:LYVE1

Note: This product is for in vitro research use only