

Lactadherin

Catalog No: #AP85349



Package Size: #AP85349-1 50ug #AP85349-2 100ug #AP85349-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Lactadherin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	MFGM, Milk fat globule-EGF factor 8, MFG-E8, SED1, Sperm surface protein SP47, MP47
Accession No.	Uniprot ID:P21956
Uniprot	P21956
GeneID	17304;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Mfge8

Note: This product is for in vitro research use only