

Matrix metalloproteinase-16

Catalog No: #AP85351



Package Size: #AP85351-1 50ug #AP85351-2 100ug #AP85351-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Matrix metalloproteinase-16
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	MMP-16,MMP-X2,Membrane-type matrix metalloproteinase 3,MT-MMP 3,Membrane-type-3 matrix metalloproteinase,MT3-MMP,MMPX2
Accession No.	Uniprot ID:P51512
Uniprot	P51512
GeneID	4325;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MMP16

Note: This product is for in vitro research use only