

Matrilysin

Catalog No: #AP85352



Package Size: #AP85352-1 50ug #AP85352-2 100ug #AP85352-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Matrilysin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Matrin, Matrix metalloproteinase-7, MMP-7, Pump-1 protease, Uterine metalloproteinase, MPSL1, PUMP1
Accession No.	Uniprot ID:P09237
Uniprot	P09237
GeneID	4316;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: MMP7

Note: This product is for in vitro research use only