

Keratin, type I cytoskeletal 17

Catalog No: #AP85353



Package Size: #AP85353-1 50ug #AP85353-2 100ug #AP85353-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Keratin, type I cytoskeletal 17
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Cytokeratin-17,CK-17,Keratin-17,K17
Accession No.	Uniprot ID:A1L595
Uniprot	A1L595
GeneID	281889;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:KRT17

Note: This product is for in vitro research use only