

Matrix metalloproteinase-14

Catalog No: #AP85355



Package Size: #AP85355-1 50ug #AP85355-2 100ug #AP85355-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Matrix metalloproteinase-14
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	MMP-14, Membrane-type matrix metalloproteinase 1, MT-MMP 1, Membrane-type-1 matrix metalloproteinase, MT-MMP, Mtmmp
Accession No.	Uniprot ID:Q10739
Uniprot	Q10739
GeneID	81707;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Mmp14

Note: This product is for in vitro research use only