

## Metallothionein-1A

Catalog No: #AP85358



Package Size: #AP85358-1 50ug #AP85358-2 100ug #AP85358-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |  |
|--------------------|--|
| Product Name       | Metallothionein-1A                                     |
| Brief Description  | Recombinant Protein                                    |
| Host Species       | E.coli   |
| Purification       | Greater than 90% by SDS-PAGE                           |
| Species Reactivity | Hu   |
| Other Names        | MT-1A, Metallothionein-IA, MT-IA, MT1S                 |
| Accession No.      | Uniprot ID:P04731                                      |
| Uniprot            | P04731   |
| GeneID             | 4489;  |
| Tag Info           | His  |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.) |

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: MT1A

Note: This product is for in vitro research use only