

Metallothionein-1E

Catalog No: #AP85359



Package Size: #AP85359-1 50ug #AP85359-2 100ug #AP85359-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Metallothionein-1E
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	MT-1E,Metallothionein-IE,MT-IE
Accession No.	Uniprot ID:P79431
Uniprot	P79431
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MT1E

Note: This product is for in vitro research use only