

MHC class I polypeptide-related sequence A

Catalog No: #AP85360



Package Size: #AP85360-1 50ug #AP85360-2 100ug #AP85360-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MHC class I polypeptide-related sequence A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	MIC-A,PERB11.1
Accession No.	Uniprot ID:Q29983
Uniprot	Q29983
GeneID	100507436;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: MICA

Note: This product is for in vitro research use only