

## Metallothionein-2

Catalog No: #AP85361



Package Size: #AP85361-1 50ug #AP85361-2 100ug #AP85361-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Metallothionein-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	MT-2, Metallothionein-II, MT-II, MT2, B28
Accession No.	Uniprot ID:Q9XST5
Uniprot	Q9XST5
GeneID	403768;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: MT2A

Note: This product is for in vitro research use only