

Mdm2-binding protein

Catalog No: #AP85362



Package Size: #AP85362-1 50ug #AP85362-2 100ug #AP85362-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mdm2-binding protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	hMTBP
Accession No.	Uniprot ID:Q96DY7
Uniprot	Q96DY7
GeneID	27085;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: MTBP

Note: This product is for in vitro research use only