

Microtubule-associated proteins 1A/1B light chain 3B

Catalog No: #AP85365



Package Size: #AP85365-1 50ug #AP85365-2 100ug #AP85365-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Microtubule-associated proteins 1A/1B light chain 3B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Autophagy-related protein LC3 B, Autophagy-related ubiquitin-like modifier LC3 B, MAP1 light chain 3-like protein 2, MAP1A/MAP1B light chain 3 B, MAP1A/MAP1B LC3 B, Microtubule-associated protein 1 light chain 3 beta, MAP1ALC3, MAP1LC3
Accession No.	Uniprot ID:O41515
Uniprot	O41515
GeneID	408001;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: MAP1LC3B

Note: This product is for in vitro research use only