

Molybdenum cofactor sulfurase

Catalog No: #AP85366



Package Size: #AP85366-1 50ug #AP85366-2 100ug #AP85366-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Molybdenum cofactor sulfurase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	MCS,Molybdenum cofactor sulfurtransferase
Accession No.	Uniprot ID:Q96EN8
Uniprot	Q96EN8
GeneID	55034;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MOCOS

Note: This product is for in vitro research use only