

Monocarboxylate transporter 2

Catalog No: #AP85367



Package Size: #AP85367-1 50ug #AP85367-2 100ug #AP85367-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Monocarboxylate transporter 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	MCT 2,Solute carrier family 16 member 7,Mct2
Accession No.	Uniprot ID:O70451
Uniprot	O70451
GeneID	20503;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Slc16a7

Note: This product is for in vitro research use only