

Nucleoside diphosphate kinase A

Catalog No: #AP85373



Package Size: #AP85373-1 50ug #AP85373-2 100ug #AP85373-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Nucleoside diphosphate kinase A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	NDK A, Metastasis inhibition factor NM23, NDPK-A, Tumor metastatic process-associated protein, nm23-M1, Nm23
Accession No.	Uniprot ID:P15532
Uniprot	P15532
GeneID	18102;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Nme1

Note: This product is for in vitro research use only