

Latent-transforming growth factor beta-binding protein

4

Catalog No: #AP85376

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP85376-1 50ug #AP85376-2 100ug #AP85376-3 1mg

Description

Product Name	Latent-transforming growth factor beta-binding protein 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	LTBP-4
Accession No.	Uniprot ID:Q8K4G1
Uniprot	Q8K4G1
GeneID	108075;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ltbp4

Note: This product is for in vitro research use only