

Lysozyme C

Catalog No: #AP85377



Package Size: #AP85377-1 50ug #AP85377-2 100ug #AP85377-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Lysozyme C
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	1,4-beta-N-acetylmuramidase C,LZM
Accession No.	Uniprot ID:P61626
Uniprot	P61626
GeneID	4069;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:LYZ

Note: This product is for in vitro research use only