

Myelin-associated glycoprotein

Catalog No: #AP85378



Package Size: #AP85378-1 50ug #AP85378-2 100ug #AP85378-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Myelin-associated glycoprotein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Siglec-4a
Accession No.	Uniprot ID:P20917
Uniprot	P20917
GeneID	17136;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Mag

Note: This product is for in vitro research use only