

Mesoderm-specific transcript protein

Catalog No: #AP85386



Package Size: #AP85386-1 50ug #AP85386-2 100ug #AP85386-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mesoderm-specific transcript protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Paternally-expressed gene 1 protein,Peg1,121a
Accession No.	Uniprot ID:Q07646
Uniprot	Q07646
GeneID	17294;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Mest

Note: This product is for in vitro research use only