

# Kin of IRRE-like protein 1

Catalog No: #AP85390



Package Size: #AP85390-1 50ug #AP85390-2 100ug #AP85390-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Kin of IRRE-like protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Kin of irregular chiasm-like protein 1,Nephrin-like protein 1,Kirrel1,Neph1
Accession No.	Uniprot ID:Q80W68
Uniprot	Q80W68
GeneID	170643;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Kirrel

Note: This product is for in vitro research use only