

## Mitochondrial carnitine/acylcarnitine carrier protein

Catalog No: #AP85392



Package Size: #AP85392-1 50ug #AP85392-2 100ug #AP85392-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |   |
|--------------------|---|
| Product Name       | Mitochondrial carnitine/acylcarnitine carrier protein                           |
| Brief Description  | Recombinant Protein   |
| Host Species       | E.coli  |
| Purification       | Greater than 90% by SDS-PAGE  |
| Species Reactivity | Rt  |
| Other Names        | Carnitine/acylcarnitine translocase,CAC,Solute carrier family 25 member 20,Cact |
| Accession No.      | Uniprot ID:P97521   |
| Uniprot            | P97521  |
| Tag Info           | His   |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.)                          |

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Slc25a20

Note: This product is for in vitro research use only