

## Protein LYRIC

Catalog No: #AP85395



Package Size: #AP85395-1 50ug #AP85395-2 100ug #AP85395-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Protein LYRIC
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Lysine-rich CEACAM1 co-isolated protein, Metadherin, Metastasis adhesion protein, Lyric
Accession No.	Uniprot ID:Q9Z1W6
Uniprot	Q9Z1W6
GeneID	170910;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Mtdh

Note: This product is for in vitro research use only