

Transcription factor Maf

Catalog No: #AP85396



Package Size: #AP85396-1 50ug #AP85396-2 100ug #AP85396-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Transcription factor Maf
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Proto-oncogene c-Maf, Transcription factor Maf-2, V-maf musculoaponeurotic fibrosarcoma oncogene homolog, Maf2
Accession No.	Uniprot ID:P54844
Uniprot	P54844
GeneID	54267;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Maf

Note: This product is for in vitro research use only