

Beta-mannosidase

Catalog No: #AP85397



Package Size: #AP85397-1 50ug #AP85397-2 100ug #AP85397-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Beta-mannosidase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Lysosomal beta A mannosidase,Mannanase,Mannase
Accession No.	Uniprot ID:Q4FZV0
Uniprot	Q4FZV0
GeneID	310864;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Manba

Note: This product is for in vitro research use only