

Mesencephalic astrocyte-derived neurotrophic factor

Catalog No: #AP85398



Package Size: #AP85398-1 50ug #AP85398-2 100ug #AP85398-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mesencephalic astrocyte-derived neurotrophic factor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Arginine-rich protein,Protein ARMET,ARMET,ARP
Accession No.	Uniprot ID:P80513
Uniprot	P80513
GeneID	540432;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MANF

Note: This product is for in vitro research use only