

Mannose-binding protein C

Catalog No: #AP85399



Package Size: #AP85399-1 50ug #AP85399-2 100ug #AP85399-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mannose-binding protein C
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	MBP-C, Mannan-binding protein, Ra-reactive factor polysaccharide-binding component p28A, RaRF p28A
Accession No.	Uniprot ID:P08661
Uniprot	P08661
GeneID	100911854;64668;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Mbl2

Note: This product is for in vitro research use only