

FGFR2 (Phospho-Ser782) Antibody

Catalog No: #12724



Package Size: #12724-1 50ul #12724-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	FGFR2 (Phospho-Ser782) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	WB
Species Reactivity	Hu Ms
Specificity	FGFR2 (Phospho-Ser782) Antibody detects endogenous levels of FGFR2 only when phosphorylated at Ser782
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human FGFR2 (Phospho-Ser782)
Target Name	FGFR2
Modification	Phospho
Other Names	FGFR2, BEK, BBDS, Bacteria-expressed kinase, CD332 antigen, CFD1, Craniofacial dysostosis 1, CEK3, Fgf bek, FGF receptor, Hydroxyaryl-protein kinase, Jackson-Weiss syndrome, JWS, FGF Receptor 2, KSAM, KGFR, Soluble FGFR4 variant 4, BFR-1, TK14, TK25, ...
Accession No.	Swiss-Prot#: P21802NCBI Gene ID: 2263
Uniprot	P21802
GeneID	2263;
Target Species	human
Calculated MW	92kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000

Note: This product is for in vitro research use only