

Monocarboxylate transporter 7

Catalog No: #AP85401



Package Size: #AP85401-1 50ug #AP85401-2 100ug #AP85401-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Monocarboxylate transporter 7
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	MCT 7, Monocarboxylate transporter 6, MCT 6, Solute carrier family 16 member 6, MCT6, MCT7
Accession No.	Uniprot ID:O15403
Uniprot	O15403
GeneID	9120;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: SLC16A6

Note: This product is for in vitro research use only