

Dual specificity mitogen-activated protein kinase kinase 6

Catalog No: #AP85402

Package Size: #AP85402-1 50ug #AP85402-2 100ug #AP85402-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Dual specificity mitogen-activated protein kinase kinase 6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	MAP kinase kinase 6,MAPK/ERK kinase 6,MEK 6
Accession No.	Uniprot ID:Q5E9X2
Uniprot	Q5E9X2
GeneID	286883;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MAP2K6

Note: This product is for in vitro research use only