

Dual specificity mitogen-activated protein kinase kinase 7

Catalog No: #AP85403

Package Size: #AP85403-1 50ug #AP85403-2 100ug #AP85403-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Dual specificity mitogen-activated protein kinase kinase 7
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	MAP kinase kinase 7, JNK-activating kinase 2, MAPK/ERK kinase 7, MEK 7, c-Jun N-terminal kinase kinase 2, JNK kinase 2, Mkk7
Accession No.	Uniprot ID:Q8CE90
Uniprot	Q8CE90
GeneID	26400;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Map2k7

Note: This product is for in vitro research use only