

# Mesoderm-specific transcript homolog protein

Catalog No: #AP85406



Package Size: #AP85406-1 50ug #AP85406-2 100ug #AP85406-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mesoderm-specific transcript homolog protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Accession No.	Uniprot ID:Q6P5P5
Uniprot	Q6P5P5
GeneID	58827;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Mest

Note: This product is for in vitro research use only