

## Sodium/calcium exchanger 1

Catalog No: #AP85409



Package Size: #AP85409-1 50ug #AP85409-2 100ug #AP85409-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Sodium/calcium exchanger 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Na(+)/Ca(2+)-exchange protein 1,Solute carrier family 8 member 1
Accession No.	Uniprot ID:P48765
Uniprot	P48765
GeneID	337925;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:SLC8A1

Note: This product is for in vitro research use only