

Myosin-1

Catalog No: #AP85411



Package Size: #AP85411-1 50ug #AP85411-2 100ug #AP85411-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Myosin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Myosin heavy chain 1,Myosin heavy chain 2x,MyHC-2x,Myosin heavy chain, skeletal muscle, adult 1
Accession No.	Uniprot ID:Q9BE40
Uniprot	Q9BE40
GeneID	281337;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MYH1

Note: This product is for in vitro research use only