

Unconventional myosin-Ib

Catalog No: #AP85413



Package Size: #AP85413-1 50ug #AP85413-2 100ug #AP85413-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Unconventional myosin-Ib
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	MYH-1c, Myosin I alpha, MMI-alpha
Accession No.	Uniprot ID:O43795
Uniprot	O43795
GeneID	4430;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: MYO1B

Note: This product is for in vitro research use only