

## Myosin-binding protein C, cardiac-type

Catalog No: #AP85414



Package Size: #AP85414-1 50ug #AP85414-2 100ug #AP85414-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Myosin-binding protein C, cardiac-type
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Cardiac MyBP-C,C-protein, cardiac muscle isoform
Accession No.	Uniprot ID:Q14896
Uniprot	Q14896
GeneID	4607;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: MYBPC3

Note: This product is for in vitro research use only