

Macrophage-expressed gene 1 protein

Catalog No: #AP85417



Package Size: #AP85417-1 50ug #AP85417-2 100ug #AP85417-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Macrophage-expressed gene 1 protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Macrophage gene 1 protein,Protein MPS1
Accession No.	Uniprot ID:A1L314
Uniprot	A1L314
GeneID	17476;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Mpeg1

Note: This product is for in vitro research use only