

MAGUK p55 subfamily member 6

Catalog No: #AP85418



Package Size: #AP85418-1 50ug #AP85418-2 100ug #AP85418-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MAGUK p55 subfamily member 6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	DIgh4 protein,P55T protein,Protein associated with Lin-7 2,DIgh4,Pals2
Accession No.	Uniprot ID:Q9JLB0
Uniprot	Q9JLB0
GeneID	56524;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Mpp6

Note: This product is for in vitro research use only