

## FGFR3 (Phospho-Tyr724) Antibody

Catalog No: #12726

Package Size: #12726-1 50ul #12726-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	FGFR3 (Phospho-Tyr724) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	FGFR3 (Phospho-Tyr724) Antibody detects endogenous levels of FGFR3 only when phosphorylated at Tyr724
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human FGFR3 (Phospho-Tyr724)
Target Name	FGFR3
Modification	Phospho
Other Names	FGFR3, CD333, CD333 antigen, ACH, FGF receptor 3, HSPGFR3EX, Hydroxyaryl-protein kinase, JTK4, FGFR-3, CEK2, Tyrosine kinase JTK4
Accession No.	Swiss-Prot#: P22607NCBI Gene ID: 2261
Uniprot	P22607
GeneID	2261;
Target Species	human
Calculated MW	95kd
Concentration	1.0mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

## Application Details

Western blotting: 1:1000

## Images



Western Blot analysis of lysates of SY-SY5Y cell, using primary antibody at 1:1000 dilution.

Note: This product is for in vitro research use only