

## Mannose-binding protein C

Catalog No: #AP85420



Package Size: #AP85420-1 50ug #AP85420-2 100ug #AP85420-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mannose-binding protein C
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	MBP-C,27 kDa mannan-binding protein monomeric subunit,pMBP-27,Mannose-binding lectin,MBL
Accession No.	Uniprot ID:Q9XSW3
Uniprot	Q9XSW3
GeneID	397230;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:MBL2

Note: This product is for in vitro research use only