

Tyrosine-protein kinase Mer

Catalog No: #AP85423



Package Size: #AP85423-1 50ug #AP85423-2 100ug #AP85423-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tyrosine-protein kinase Mer
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Proto-oncogene c-Mer, Receptor tyrosine kinase MerTK, Mer
Accession No.	Uniprot ID:Q60805
Uniprot	Q60805
GeneID	17289;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: MERTK

Note: This product is for in vitro research use only