

S-adenosylmethionine synthase isoform type-2

Catalog No: #AP85429



Package Size: #AP85429-1 50ug #AP85429-2 100ug #AP85429-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	S-adenosylmethionine synthase isoform type-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	AdoMet synthase 2, Methionine adenosyltransferase 2, MAT 2
Accession No.	Uniprot ID:Q3THS6
Uniprot	Q3THS6
GeneID	232087;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Mat2a

Note: This product is for in vitro research use only