

## S-adenosylmethionine synthase isoform type-1

Catalog No: #AP85430



Package Size: #AP85430-1 50ug #AP85430-2 100ug #AP85430-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	S-adenosylmethionine synthase isoform type-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	AdoMet synthase 1,Methionine adenosyltransferase 1,MAT 1
Accession No.	Uniprot ID:Q2KJC6
Uniprot	Q2KJC6
GeneID	541078;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:MAT1A

Note: This product is for in vitro research use only