

Methylated-DNA--protein-cysteine methyltransferase

Catalog No: #AP85431



Package Size: #AP85431-1 50ug #AP85431-2 100ug #AP85431-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Methylated-DNA--protein-cysteine methyltransferase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	6-O-methylguanine-DNA methyltransferase, MGMT, O-6-methylguanine-DNA-alkyltransferase
Accession No.	Uniprot ID:P26187
Uniprot	P26187
GeneID	17314;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Mgmt

Note: This product is for in vitro research use only