

CDK-activating kinase assembly factor MAT1

Catalog No: #AP85435



Package Size: #AP85435-1 50ug #AP85435-2 100ug #AP85435-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | CDK-activating kinase assembly factor MAT1 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Ms |
| Other Names | CDK7/cyclin-H assembly factor, Menage a trois, RING finger protein MAT1, p35, p36, Mat1 |
| Accession No. | Uniprot ID:P51949 |
| Uniprot | P51949 |
| GeneID | 17420; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Mnat1

Note: This product is for in vitro research use only