

## Canalicular multispecific organic anion transporter 2

Catalog No: #AP85436



Package Size: #AP85436-1 50ug #AP85436-2 100ug #AP85436-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Canalicular multispecific organic anion transporter 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	ATP-binding cassette sub-family C member 3,Multi-specific organic anion transporter D,MOAT-D,Multidrug resistance-associated protein 3,CMOAT2,MLP2,MRP3
Accession No.	Uniprot ID:O15438
Uniprot	O15438
GeneID	8714;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:ABCC3

Note: This product is for in vitro research use only