

Neogenin

Catalog No: #AP85443



Package Size: #AP85443-1 50ug #AP85443-2 100ug #AP85443-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Neogenin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Immunoglobulin superfamily DCC subclass member 2,IGDCC2,NGN
Accession No.	Uniprot ID:Q92859
Uniprot	Q92859
GeneID	4756;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: NEO1

Note: This product is for in vitro research use only